

Rat IgG1, kappa Isotype Control Antibody (RTK2071)

Summary

Catalog No.	BGT-ABS-451
Host species	Rat
Conjugation	Unconjugated
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95%
Clonality	Monoclonal
Isotype	IgG1, kappa
Applications	FCM,Control
Purification	Protein A or G purified.
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.Store at +4°C short term (1-2 weeks).Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only.