

7-methylguanosine (m7G) Monoclonal Antibody

Overview

Catalog No.	BGT-ANT-03172
Host species	Mouse
Form	Liquid
Storage buffer	0.01M PBS pH 7.4
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a, kappa
Applications	ELISA, DB, IP, IHC
Target	7-methylguanosine, m7G
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Clone ID	SAA2426
Note	For research use only