

Human CXCL4/PF4 Antibody, Clone 5B9

Overview

Catalog No.	BGT-ANT-03185	
Host species	Mouse	
Species reactivity	Human	
Form	Liquid	
Storage buffer	0.01M PBS, pH 7.4.	
Purity	>95% as determined by SDS-PAGE.	
Clonality	Monoclonal	
Isotype	IgG2a, kappa	
Applications	ELISA	
Target	C-X-C motif chemokine 4, PF-4, Oncostatin-A, Endothelial cell growth inhibitor, CXCL4, Iroplact, Platelet factor 4, PF4, SCYB4	
Purification	Protein A/G purified from cell culture supernatant.	
Endotoxin level	Please contact with the lab for this information.	
Accession	P02776	
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.	
Clone ID	5B9	
Product Usage Information	Application	Dilution
	ELISA	1:2000-1:20000
Note	For research use only.	
RRID	AB_3729716	