


AQP9 Polyclonal antibody

Summary

Catalog No.	BGT-ANT-10267
Alternative Names	AQP-9,SSC1,Small solute channel 1,Aquaporin-9,AQP9,Aquaglyceroporin-9
Host species	Rabbit
Tested applications	ELISA: 1:4000-1:8000, IHC: 1:50-1:100, WB: 1:1000-1:4000
Species reactivity	Human
Immunogen	E. coli - derived recombinant Human AQP9 (Tyr125-Leu158).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	WB,IHC,ELISA
Purification	Antigen affinity purification.
Accession	O43315
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 4°C for frequent use.Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only.

Data Image



WESTERN BLOT

Recombinant Protein lysates were subjected to SDS PAGE followed by western blot with AQP9 antibody (BGT-ANT-10267) at 1 $\mu\text{g/ml}$.

Lane 1: Recombinant protein

Second Ab: Goat Anti-Rabbit IgG H&L Polyclonal antibody, HRP at 0.1 $\mu\text{g/mL}$.

Predict MW: 30 kDa

Observed MW: 30 kDa