

ASFV p72/B646L Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-10324
Alternative Names	Hexon protein p72, Major capsid protein, MCP, p72, p73, Ba71V-081, B646L, African swine fever virus (ASFV)
Host species	Rabbit
Species reactivity	African swine fever virus (ASFV)
Immunogen	E. coli - derived recombinant ASFV p72/B646L (Ile423-Val641).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, WB
Purification	Antigen affinity purification.
Accession	P22776
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 4°C for frequent use.Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only.