

BAMBI Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-10405
Alternative Names	BAMBI, BMP and activin membrane-bound inhibitor homolog, Non-metastatic gene A protein, NMA, Putative transmembrane protein NMA
Host species	Rabbit
Tested applications	ELISA: 1:4000-1:8000, IHC: 1:50-1:100, WB: 1:1000-1:4000
Species reactivity	Human, Mouse, Rat
Immunogen	E. coli - derived recombinant Human BAMBI (Arg29-Asp130).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA,IHC,WB
Target	BAMBI
Purification	Antigen affinity purification.
Accession	Q13145
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 4°C for frequent use.Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only.