

BVDV E2 glycoprotein Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-10495
Alternative Names	Genome polyprotein, N-terminal protease, N-pro, 3.4.22.-, Autoprotease p20, Capsid protein C, E(rns) glycoprotein, gp44/48, Envelope glycoprotein E1, gp33, Envelope glycoprotein E2, gp55, Viroporin p7, Non-structural protein 2-3, Cysteine protease NS2, 3.4.22.-, Non-structural protein 2, Serine protease NS3, 3.4.21.113, 3.6.1.15, 3.6.4.13, Non-structural protein 3, Non-structural protein 4A, NS4A, Non-structural protein 4B, NS4B, Non-structural protein 5A, NS5A, RNA-directed RNA polymerase, 2.7.7.48, NS5B, Bovine viral diarrhea virus (BVDV)
Host species	Rabbit
Species reactivity	Bovine viral diarrhea virus (BVDV) (Mucosal disease virus)
Immunogen	E. coli - derived recombinant BVDV (Mucosal disease virus) E2 glycoprotein (Cys1-Asp218).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, WB
Purification	Antigen affinity purification.
Accession	AGC12054.1

Stability and Storage

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80 °C for twelve months from the date of receipt.

Note

For research use only.