

BVDV Envelope glycoprotein E2 Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-10496
Alternative Names	Genome polyprotein, Envelope glycoprotein E2, gp55, Bovine viral diarrhea virus (isolate NADL) (BVDV) (Mucosal disease virus)
Host species	Rabbit
Species reactivity	Bovine viral diarrhea virus (isolate NADL) (BVDV) (Mucosal disease virus)
Immunogen	E. coli - derived recombinant BVDV Envelope glycoprotein E2 (His693-Ser851).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, WB
Purification	Antigen affinity purification.
Accession	P19711
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 4°C for frequent use.Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only.