

CD16/CD32 Monoclonal Antibody

Overview

Catalog No.	BGT-ANT-10816
Host species	Rat
Species reactivity	Mouse
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2b, kappa
Applications	FCM, FuncS, IF, IP
Target	Low affinity immunoglobulin gamma Fc region receptor III, IgG Fc receptor III, Fc-gamma RIII, FcRIII, CD16, Low affinity immunoglobulin gamma Fc region receptor II, Fc gamma receptor IIB, Fc-gamma RII, Fc-gamma-RIIB, FcRII, IgG Fc receptor II beta, Lymphocyte antigen 17, Ly-17, CD32, Fcgr2, Fcgr2b, Ly-17
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	P08508, P08101
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Clone ID	2.4G2
Note	For research use only.