

CD16 Monoclonal Antibody

Summary

Catalog No.	BGT-ANT-10818
Alternative Names	CD16, CD16a, CD16A, FCGR3A, Fc-gamma RIII-alpha, FcRIIIa, IgG Fc receptor III-2, Low affinity immunoglobulin gamma Fc region receptor III-A, IGFR3, CD16a antigen, Fc-gamma RIII, FCGR3, Fc-gamma RIIIa, FCG3, FcRIII, FcR-10, Low affinity immunoglobulin gamma Fc region receptor III-B, Fc-gamma RIII-beta, CD16-I, Fc-gamma RIII, Fc-gamma RIIIb, FcRIII, FcRIIIb, FcR-10, IgG Fc receptor III-1, CD16b, FCGR3B, CD16B, FCG3, FCGR3, IGFR3
Host species	Mouse
Conjugation	Unconjugated
Species reactivity	Human
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95%
Clonality	Monoclonal
Isotype	IgG1, kappa
Applications	Cell Isolation, CyTOF, Depletion, ELISA, FCM, Functional Assay, ICC, IF, IHC, IP, Neutralization and Blocking
Purification	Protein A or G purified.
Endotoxin level	Please contact with the lab for this information.

Accession P08637 & O75015

Stability and Storage Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.

Note For research use only.