


CD200R1 / OX2R Monoclonal Antibody

Summary

Catalog No.	BGT-ANT-10905
Alternative Names	MOX2R, Cell surface glycoprotein OX2 receptor 1, OX2R, Cell surface glycoprotein CD200 receptor 1, CD200R, CD200 cell surface glycoprotein receptor, CRTR2, CD200R1
Host species	Mouse
Conjugation	Unconjugated
Species reactivity	Human
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	2.96 mg/ml
Purity	>95%
Clonality	Monoclonal
Isotype	IgG1, kappa
Applications	ELISA, FCM
Purification	Protein A or G purified.
Endotoxin level	Please contact with the lab for this information.
Accession	Q8TD46
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.Store at +4°C short term (1-2 weeks).Store at -20 °C 12 months. Store at -80°C long term.

Note For research use only.

Data Image



Detects CD200R1/OX2R in indirect ELISAs.

Bioactivity