

CD200R1 / OX2R Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-10906
Alternative Names	Cell surface glycoprotein CD200 receptor 1 (CD200 cell surface glycoprotein receptor) (Cell surface glycoprotein OX2 receptor 1)
Host species	Rabbit
Species reactivity	Mouse
Immunogen	E. coli - derived recombinant Mouse CD200R1/OX2R (Gln42-Ser229).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Purification	Antigen affinity purification.
Accession	Q9ES57
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 4°C for frequent use.Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only.