

CD25 / IL2RA Monoclonal Antibody(PC61/PC61.5.3)

Summary

Catalog No.	BGT-ANT-10998
Alternative Names	Interleukin-2 receptor subunit alpha, IL-2 receptor subunit alpha, IL-2-RA, IL-2R subunit alpha, IL2-RA, p55, CD25, Il2ra, Il2r
Host species	Rat
Conjugation	Unconjugated
Species reactivity	Mouse
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	2.01 mg/ml
Purity	>95%
Clonality	Monoclonal
Isotype	IgG1, lambda
Applications	FCM
Purification	Protein A or G purified.
Endotoxin level	Please contact with the lab for this information.
Accession	P01590
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.Store at +4°C short term (1-2 weeks).Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only.