

CD279 / PD1 / PD1 Monoclonal Antibody

Summary

Catalog No.	BGT-ANT-11044
Alternative Names	Programmed cell death protein 1, Protein PD-1, hPD-1, PD1, PDCD1, CD279
Host species	Rat
Conjugation	Unconjugated
Species reactivity	Mouse
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a, kappa
Applications	FCM, IF, IHC, Neutralization, WB
Purification	Protein A or G purified.
Endotoxin level	Please contact with the lab for this information.
Accession	Q02242
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.Store at +4°C short term (1-2 weeks).Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only.