

CD354/TREM1 Polyclonal Antibody

Overview

Catalog No.	BGT-ANT-11166
Host species	Rabbit
Species reactivity	Mouse
Immunogen	E. coli - derived recombinant Mouse CD354/TREM1 (Glu26-Lys135).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Target	TREM1, Triggering receptor expressed on monocytes 1, Triggering receptor expressed on myeloid cells 1, TREM-1, CD354
Purification	Purified by antigen affinity column.
Accession	Q9JKE2
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

	Application	Dilution
Product Usage Information	ELISA	1:5000-1:20000
	IHC	1:50-1:500
	WB	1:500-1:2000
Note	For research use only.	