

CLDN18/Claudin-18 Polyclonal Antibody

Overview

Catalog No.	BGT-ANT-11535	
Host species	Rabbit	
Species reactivity	Human, Mouse, Rat	
Immunogen	E. coli - derived recombinant Human CLDN18/Claudin-18 (Thr25-Gly65).	
Form	Liquid	
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.	
Clonality	Polyclonal	
Isotype	IgG	
Applications	ELISA, IHC, WB	
Target	Claudin-18, CLDN18	
Purification	Purified by antigen affinity column.	
Accession	P56856	
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.	
Product Usage Information	Application	Dilution
	ELISA	1:5000-1:20000
	IHC	1:50-1:500
	WB	1:500-1:2000
Note	For research use only.	