

CLDN3 / 4 Monoclonal Antibody

Summary

Catalog No.	BGT-ANT-11537
Alternative Names	Claudin-3, Rat ventral prostate.1 protein homolog, CPE-R 2, C7orf1, Clostridium perfringens enterotoxin receptor 2, CPE-receptor 2, CPETR2, hRVP1, CLDN3, CPER, CLDN4, Williams-Beuren syndrome chromosomal region 8 protein, CPE-R, CPE-receptor, WBSCR8, Claudin-4, Clostridium perfringens enterotoxin receptor, CPETR1
Host species	Mouse
Conjugation	Unconjugated
Species reactivity	Human
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a, kappa
Applications	ELISA, FCM
Purification	Protein A or G purified.
Endotoxin level	Please contact with the lab for this information.
Accession	O15551, O14493

Stability and Storage

Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.

Note

For research use only.