

CXCL13 / BCA-1 / BLC Antibody

Summary

Catalog No.	BGT-ANT-11763
Clone ID	SAA0468
Host species	Mouse
Conjugation	Unconjugated
Species reactivity	Human
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a, lambda
Applications	ELISA, FCM, Neutralization
Target	SCYB13, BLC, CXCL13, B cell-attracting chemokine 1, B lymphocyte chemoattractant, Small-inducible cytokine B13, CXC chemokine BLC, Angie, C-X-C motif chemokine 13, BCA-1, BCA1
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	O43927

Stability and Storage

Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Note

For research use only.