


EBF1 Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-11956
Alternative Names	EBF,EBF1,OE-1,COE1,Early B-cell factor,Transcription factor COE1,O/E-1,
Host species	Rabbit
Tested applications	WB:1:2000~1:8000, IHC:1:50~100 , Elisa:1:4000~1:8000
Species reactivity	Human
Immunogen	E. coli - derived recombinant Human EBF1 (Asp179-Ser451).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	WB,IHC,ELISA
Purification	Antigen affinity purification.
Accession	Q9UH73
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 4°C for frequent use.Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only.

Data Image



Western blot

Various lysates were subjected to SDS PAGE followed by western blot with EBF1 antibody (BGT-ANT-11956) at 0.42µg/ml.

Lane 1: Raji cell lysate

Lane 2: THP-1 cell lysate

Lane 3: K562 cell lysate

Second Ab: Goat Anti-Rabbit IgG H&L Polyclonal antibody, HRP at 0.1 µg/mL.

Predict MW: 64 kDa

Observed MW: 64 kDa