


EFNA1 Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-11982
Host species	Rabbit
Tested applications	ELISA: 1:4000-1:8000, IHC: 1:50-1:100, WB: 1:1000-1:4000
Species reactivity	Human, Mouse, Rat
Immunogen	E. coli - derived recombinant Human EFNA1 (Asp19-Ser182).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Target	EPLG1, EFNA1, LERK1, Immediate early response protein B61, EPH-related receptor tyrosine kinase ligand 1, TNFAIP4, TNF alpha-induced protein 4, LERK-1, Ephrin-A1, Tumor necrosis factor alpha-induced protein 4
Purification	Purified by antigen affinity column.
Accession	P20827
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only.

Data Image



Western Blot

Recombinant Protein lysates were subjected to SDS PAGE followed by western blot with EFNA1 antibody (BGT-ANT-11982) at 1 μ g/ml.

Lane 1: Recombinant protein

Second Ab: Goat Anti-Rabbit IgG H&L Polyclonal antibody, HRP at 0.1 μ g/mL.

Predict MW: 22 kDa

Observed MW: 22 kDa