

Human FGF23 Monoclonal Antibody

Overview

Catalog No.	BGT-ANT-12185
Host species	Mouse
Species reactivity	Human
Form	Liquid
Storage buffer	0.01M PBS pH 7.4
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a, kappa
Applications	ELISA, FCM
Target	Phosphatonin, FGF-23, HYPF, FGF23, Fibroblast growth factor 23, Tumor-derived hypophosphatemia-inducing factor
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	Q9GZV9
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Clone ID	SAA0441
Note	For research use only.