

GRIN1 / NMDAR1 Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-12417
Alternative Names	N-methyl-D-aspartate receptor subunit NR1, GluN1, Glutamate [NMDA] receptor subunit zeta-1, Glutamate receptor ionotropic, NMDA 1, NMD-R1, GRIN1, NMDAR1
Host species	Rabbit
Species reactivity	Human
Immunogen	E. coli - derived recombinant Human GRIN1/NMDAR1 (Asp23-His371).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Purification	Antigen affinity purification.
Accession	Q05586
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 4°C for frequent use.Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only.