

GSDME Polyclonal Antibody

Overview

Catalog No.	BGT-ANT-12429	
Host species	Rabbit	
Species reactivity	Human	
Immunogen	E. coli - derived recombinant Human GSDME (Met1-Cys235).	
Form	Liquid	
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.	
Clonality	Polyclonal	
Isotype	IgG	
Applications	ELISA, IHC, WB	
Target	ICERE-1, ICERE1, GSDME-CT, GSDME-NT, DFNA5, Inversely correlated with estrogen receptor expression 1, Gasdermin-E, GSDME, Non-syndromic hearing impairment protein 5	
Purification	Purified by antigen affinity column.	
Accession	O60443	
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.	
Product Usage Information	Application	Dilution
	ELISA	1:5000-1:20000
	IHC	1:50-1:500
	WB	1:500-1:2000
Note	For research use only.	