

H2-Ab1 / MHCII Monoclonal Antibody

Summary

Catalog No.	BGT-ANT-12463
Alternative Names	H2-Ab1, H-2 class II histocompatibility antigen, A beta chain, H2-iabeta
Host species	Mouse
Conjugation	Unconjugated
Species reactivity	Mouse
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95%
Clonality	Monoclonal
Isotype	VHH, mFc
Applications	FCM,IF
Purification	Protein A or G purified.
Endotoxin level	Please contact with the lab for this information.
Accession	P14483
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.Store at +4°C short term (1-2 weeks).Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only.