


HSPB7 Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-12633
Host species	Rabbit
Tested applications	ELISA: 1:4000-1:8000, IHC: 1:50-1:100, WB: 1:1000-1:4000
Species reactivity	Human, Mouse, Rat
Immunogen	E. coli - derived recombinant Human HSPB7 (Met1-Ile170).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Target	HSPB7,CVHSP,Cardiovascular heat shock protein,cvHsp,HspB7,Heat shock protein beta-7
Purification	Purified by antigen affinity column.
Accession	Q9UBY9
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only.

Data Image



WESTERN BLOT

Recombinant Protein lysates were subjected to SDS PAGE followed by western blot with HSPB7 antibody (BGT-ANT-12633) at 1 µg/ml.

Lane 1: Recombinant protein

Second Ab: Goat Anti-Rabbit IgG H&L Polyclonal antibody, HRP at 0.1 µg/mL.

Predict MW: 21 kDa

Observed MW: 21 kDa