


MSLN Monoclonal Antibody

Summary

Catalog No.	BGT-ANT-13309
Alternative Names	CAK1 antigen, Mesothelin, MSLN, Pre-pro-megakaryocyte-potentiating factor, MPF
Host species	Mouse
Conjugation	Unconjugated
Species reactivity	Human
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95%
Clonality	Monoclonal
Isotype	IgG2a, lambda
Applications	FCM, IHC
Purification	Protein A or G purified.
Endotoxin level	Please contact with the lab for this information.
Accession	Q13421
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.Store at +4°C short term (1-2 weeks).Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only.

Data Image



Flow-cytometry using anti-human MSLN antibody. HeLa cells were stained with an irrelevant antibody (Blue Histogram) or an anti-human MSLN antibody monoclonal antibody (Catalog # BGT-ANT-13309, Green Histogram) at a concentration of 5 $\mu\text{g/ml}$ for 30 mins at RT. After washing, bound antibody was detected using a FITC conjugated goat anti-mouse antibody and cells analysed on a NovoCyt Flow Cytometer.

Flow-cytometry