

MUL1 Polyclonal Antibody


Summary

Catalog No.	BGT-ANT-13340
Host species	Rabbit
Tested applications	ELISA: 1:4000-1:8000, IHC: 1:50-1:100, WB: 1:1000-1:4000
Species reactivity	Human
Immunogen	E. coli - derived recombinant Human MUL1 (Lys44-Gln224).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Target	Mitochondrial-anchored protein ligase, E3 ubiquitin-protein ligase MUL1, Protein Hades, C1orf166, RNF218, Mitochondrial ubiquitin ligase activator of NFKB 1, Growth inhibition and death E3 ligase, MUL1, RING finger protein 218, MAPL, RING-type E3 ubiquitin transferase NFKB 1, Putative NF-kappa-B-activating protein 266, GIDE, MULAN, E3 SUMO-protein ligase MUL1
Purification	Purified by antigen affinity column.
Accession	Q969V5
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note

For research use only.

Data Image



WESTERN BLOT

Various lysates were subjected to SDS PAGE followed by western blot with MUL1 antibody (BGT-ANT-13340) at 1 μ g/ml.

Lane 1: Mouse skeletal muscle lysate

Lane 2: Rat skeletal muscle lysate

Second Ab: Goat Anti-Rabbit IgG H&L Polyclonal antibody, HRP at 0.1 μ g/mL.

Predict MW: 40 kDa

Observed MW: 40 kDa