

Mycoplasma pneumoniae mgpA/Adhesin P1

Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-13350
Host species	Rabbit
Tested applications	ELISA: 1:4000-1:8000, WB: 1:2000-1:8000
Species reactivity	Mycoplasma pneumoniae (strain ATCC 29342 / M129 / Subtype 1) (Mycoplasma pneumoniae)
Immunogen	E. coli - derived recombinant Mycoplasma pneumoniae mgpA/Adhesin P1 (Ala1160-Asp1521).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Target	Adhesin P1, Attachment protein, Cytadhesin P1, mgpA, MPN_141, MP013
Purification	Purified by antigen affinity column.
Accession	P11311
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only.