

Pirb / Lilrb3 Polyclonal Antibody

Overview

Catalog No.	BGT-ANT-13776
Host species	Rabbit
Species reactivity	Mouse
Immunogen	E. coli - derived recombinant Mouse Pirb/Lilrb3 (Ile30-Gly316).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Target	Paired immunoglobulin-like receptor B, Pirb, Lilrb3
Purification	Purified by antigen affinity column.
Accession	Q8K4V6
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only.