

PMFBP1 Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-13822
Alternative Names	Polyamine-modulated factor 1-binding protein 1, PMFBP1, PMF-1-binding protein
Host species	Rabbit
Tested Applications	ELISA: 1:4000~1:8000, IHC: 1:50~100, WB: 1:2000~1:4000
Species reactivity	Human
Immunogen	E. coli - derived recombinant Human PMFBP1 (Met1-Asn150).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Purification	Antigen affinity purification.
Accession	Q8TBY8
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 4°C for frequent use.Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only.