


PPBP / CXCL7 / NAP-2 Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-13872
Host species	Rabbit
Tested applications	ELISA: 1:4000-1:8000, IHC: 1:50-1:100, WB: 1:1000-1:4000
Species reactivity	Human
Immunogen	E. coli - derived recombinant Human PPBP/CXCL7/NAP-2 (Ala59-Asp128).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Target	NAP-2,NAP-2(1-63),NAP-2(73),PPBP,CXCL7,PBP,Small-inducible cytokine B7,Macrophage-derived growth factor,LA-PF4,Leukocyte-derived growth factor,CTAP-III,CTAP-III(1-81),Low-affinity platelet factor IV,TGB1,Beta-TG,C-X-C motif chemokine 7,MDGF,Platelet basic protein,LDGF,NAP-2(74),THBGB1,NAP-2(1-66),CTAP3,SCYB7
Purification	Purified by antigen affinity column.
Accession	P02775
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only.

Data Image



WESTERN BLOT

Recombinant Protein lysates were subjected to SDS PAGE followed by western blot with PPBP/CXCL7/NAP-2 antibody (BGT-ANT-13872) at 1 $\mu\text{g/ml}$.

Lane 1: Recombinant protein

Second Ab: Goat Anti-Rabbit IgG H&L Polyclonal antibody, HRP at 0.1 $\mu\text{g/mL}$.

Predict MW: 10 kDa

Observed MW: 10 kDa