

PRDX5 Polyclonal antibody

Summary

Catalog No.	BGT-ANT-13902
Alternative Names	Peroxisomal antioxidant enzyme,Antioxidant enzyme B166,Thioredoxin peroxidase PMP20,AOEB166,TPx type VI,Peroxiredoxin-5,mitochondrial,ACR1,Liver tissue 2D-page spot 71B,PRDX5,PLP,Thioredoxin-dependent peroxiredoxin 5,Alu corepressor 1,Peroxiredoxin V,Prx-V,
Host species	Rabbit
Tested applications	WB:1:2000~1:8000, IHC:1:50~100 , Elisa:1:4000~1:8000
Species reactivity	Human
Immunogen	E. coli - derived recombinant Human PRDX5 (Met1-Asp162).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	WB,IHC,ELISA
Purification	Antigen affinity purification.
Accession	P30044
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 4°C for frequent use.Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only.