


Rabies virus / RABV G / Glycoprotein Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-14054
Host species	Rabbit
Tested applications	ELISA: 1:4000-1:8000, WB: 1:1000-1:4000
Species reactivity	Rabies virus (strain CVS-11) (RABV)
Immunogen	E. coli - derived recombinant Rabies virus/RABV G/Glyco(Leu30-Tyr459).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Target	Glycoprotein, G
Purification	Purified by antigen affinity column.
Accession	O92284
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only.

Data Image



WESTERN BLOT

Recombinant Protein lysates were subjected to SDS PAGE followed by western blot with G/Glycoprotein antibody (BGT-ANT-14054) at 1 $\mu\text{g/ml}$.

Lane 1: Recombinant protein

Second Ab: Goat Anti-Rabbit IgG H&L Polyclonal antibody, HRP at 0.1 $\mu\text{g/mL}$.

Predict MW: 51 kDa

Observed MW: 51 kDa