

SLC20A1 Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-14420
Alternative Names	GLVR1, Gibbon ape leukemia virus receptor 1, GLVR-1, PIT1, Phosphate transporter 1, Leukemia virus receptor 1 homolog, PiT-1, Solute carrier family 20 member 1, Sodium-dependent phosphate transporter 1, SLC20A1
Host species	Rabbit
Species reactivity	Human
Immunogen	E. coli - derived recombinant Human SLC20A1 (Cys251-Glu510).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Purification	Antigen affinity purification.
Accession	Q8WUM9
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 4°C for frequent use.Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only.