

SMOC2 Polyclonal Antibody

Overview

Catalog No.	BGT-ANT-14460	
Host species	Rabbit	
Species reactivity	Mouse	
Immunogen	E. coli - derived recombinant Mouse SMOC2 (Gly306-Gly414).	
Form	Liquid	
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.	
Clonality	Polyclonal	
Isotype	IgG	
Applications	ELISA, IHC, WB	
Target	SMAP-2, Smooth muscle-associated protein 2, SMOC-2, SMAP2, SMOC2, Secreted modular calcium-binding protein 2, SPARC-related modular calcium-binding protein 2	
Purification	Purified by antigen affinity column.	
Accession	Q8CD91	
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.	
Product Usage Information	Application	Dilution
	ELISA	1:5000-1:20000
	IHC	1:50-1:500
	WB	1:500-1:2000
Note	For research use only.	