


ST6GALNAC6 Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-14556
Host species	Rabbit
Tested applications	ELISA: 1:4000-1:8000, IHC: 1:50-1:100, WB: 1:1000-1:4000
Species reactivity	Human, Mouse, Rat
Immunogen	E. coli - derived recombinant Human ST6GALNAC6 (Val67-Thr333).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Target	hST6GalNAc VI,SIAT7-F,Sialyltransferase 7F,ST6GalNAcVI,SIAT7F,Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 6,GalNAc alpha-2,6-sialyltransferase VI,ST6GalNAc VI,ST6GALNAC6
Purification	Purified by antigen affinity column.
Accession	Q969X2
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only.

Data Image



Western blot

Various lysates were subjected to SDS PAGE followed by western blot with CA19-9 antibody (BGT-ANT-14556) at 0.64ug/ml.

Lane 1: LOVO cell lysate

Second Ab: Goat Anti-Rabbit IgG H&L Polyclonal antibody, HRP at 0.1 µg/mL.

Predict MW: 34 kDa

Observed MW: 36 kDa