

Mouse ZNRF3 Polyclonal Antibody

Summary

Catalog No.	BGT-ANT-15126
Host species	Rabbit
Tested applications	ELISA: 1:4000-1:8000, IHC: 1:50-1:100, WB: 1:1000-1:4000
Species reactivity	Mouse
Immunogen	E. coli - derived recombinant Mouse ZNRF3 (Lys53-Pro207).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, IHC, WB
Target	KIAA1133, RNF203, RING-type E3 ubiquitin transferase ZNRF3, ZNRF3, RING finger protein 203, Zinc/RING finger protein 3, E3 ubiquitin-protein ligase ZNRF3
Purification	Purified by antigen affinity column.
Accession	Q5SSZ7
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only.