

4R-HNE/4-HNE Antibody

Overview

Catalog No.	BGT-ANT-35060	
Host species	Mouse	
Species reactivity	General	
Form	Liquid	
Storage buffer	0.01M PBS, pH 7.4.	
Purity	>95% as determined by SDS-PAGE.	
Clonality	Monoclonal	
Isotype	IgG2a	
Applications	ELISA, WB	
Target	HNE, 4-HNE, 4-hydroxynon-2-enal, 4-Hydroxy-2-nonenal, 4-Hydroxynonenal	
Purification	Protein A/G purified from cell culture supernatant.	
Endotoxin level	Please contact with the lab for this information.	
Accession	ChEBI: 32585, 29343-52-0	
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.	
Clone ID	R310	
Product Usage Information	Application	Dilution
	ELISA	1:2000-1:20000
	WB	1:500-1:2000
Note	For research use only.	