

ASFV B117L Antibody

Overview

Catalog No.	BGT-ANT-35102
Host species	Mouse
Species reactivity	African swine fever virus
Immunogen	Recombinant ASFV B117L protein.
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA, WB
Target	hypothetical protein IM014_gp118, B117L CDS
Purification	Protein A/G purified from cell culture supernatant.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Clone ID	1A026
Note	For research use only.