

Coxsackievirus A6/CVA6 VP1/Capsid protein VP1 Antibody

Summary

Catalog No.	BGT-ANT-35259
Clone ID	SAA0354
Host species	Human
Species reactivity	Coxsackievirus A6
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG1, kappa
Applications	ELISA, SPR
Target	Genome polyprotein, Capsid protein VP1, P1D, Virion protein 1
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	A0A0K2BNC7
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.