


HPV16 E6/Protein E6 Antibody

Summary

Catalog No.	BGT-ANT-35508
Clone ID	Nb9
Host species	Mouse
Species reactivity	Human papillomavirus type 16
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	VHH-mFC
Applications	ELISA, Neutralization, WB
Target	E6, Protein E6
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Expression system	Mammalian Cells
Accession	P03126
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Data Image



SDS-PAGE

SDS PAGE for HPV16 E6/Protein E6 Antibody