

Human oxLDL Antibody

Summary

Catalog No.	BGT-ANT-36216
Clone ID	SAA0551
Host species	Human
Species reactivity	Human
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA
Target	oxidised LDL, oxidised form of low-density lipoprotein (LDL), oxidised low-density lipoproteins, oxidized low-density lipoprotein, oxidized low-density lipoproteins, OxLDL
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.