

KCC2 Antibody

Summary

Catalog No.	BGT-ANT-36568
Alternative Names	CAMKG, CAMK-II, CAMK, CAMK2G, CaM kinase II subunit gamma, CaMK-II subunit gamma, Calcium/calmodulin-dependent protein kinase type II subunit gamma
Host species	Mouse
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a, kappa
Applications	IHC, WB
Species reactivity	Human, Mouse, Rat
Purification	Protein A or G purified from cell culture supernatant.
Endotoxin level	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.Store at +4°C short term (1-2 weeks).Store at -20 °C 12 months. Store at -80°C long term.
Stability and Storage	For research use only.
Note	Please contact with the lab for this information.