

Staphylococcal protein A/SpA Antibody

Summary

Catalog No.	BGT-ANT-37037
Clone ID	G3
Host species	Human
Species reactivity	Staphylococcus aureus
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG3, Kappa
Applications	ELISA
Target	Immunoglobulin G-binding protein A, IgG-binding protein A, Staphylococcal protein A, SpA
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	P02976
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.