

Suid herpesvirus 1/SuHV-1 UL42/PAP Antibody

Summary

Catalog No.	BGT-ANT-37051
Clone ID	2D5
Host species	Mouse
Species reactivity	Suid herpesvirus 1 (strain Kaplan) (SuHV-1) (Pseudorabies virus (strain Kaplan)), Suid alphaherpesvirus 1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2, kappa
Applications	ELISA, IF, IP, WB
Target	UL42, DNA polymerase processivity factor, Polymerase accessory protein, PAP
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	P36702, A0A0F6QNH2


Stability and Storage

Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Note

For research use only.

Data Image



SDS PAGE for Suid herpesvirus 1/SuHV-1 UL42/PAP Antibody

SDS-PAGE