


Z-DNA Antibody

Summary

Catalog No.	BGT-ANT-37165
Clone ID	Z22
Host species	Mouse
Specificity	Z22 binds both Z-DNA, but not B-DNA or ssDNA (single-stranded DNA).
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.89 mg/ml
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a, kappa
Applications	DB, ELISA, IF
Target	DNA, Z-DNA
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Data Image



SDS PAGE for Z-DNA Antibody (Z22)

SDS-PAGE