

α -Gal Epitope (Gal α 1-3Gal β 1-4GlcNAc-R) Antibody

Summary

Catalog No.	BGT-ANT-37177
Clone ID	B1H12
Host species	Mouse
Species reactivity	Mouse, Rat, Porcine
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA, FCM, FuncS, IHC-P, WB
Target	α -Galactose 1,3, α -Gal Epitope (Gal α 1-3Gal β 1-4GlcNAc-R)
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.