

EBV/HHV4 gB/BALF4 Antibody, InVivo Block/Neutralizing Antibody


Summary

Catalog No.	BGT-ANT-37212
Clone ID	Iv0133
Host species	Human
Species reactivity	Epstein-Barr virus (strain GD1) (HHV-4) (Human herpesvirus 4)
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.2 mg/ml
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG1, kappa
Applications	Neutralization, WB
Target	Envelope glycoprotein B, gB, BALF4, gB, EBVaGC1_061
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	R4R670
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Note

For research use only. Not suitable for clinical or therapeutic use.

Data Image



SDS PAGE for EBV/HHV4 gB/BALF4 Antibody

SDS-PAGE