

Human FGF21, InVivo Block/Neutralizing Antibody

Overview

Catalog No.	BGT-ANT-37337
Host species	Human
Species reactivity	Human
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG
Applications	ELISA, Neutralization
Target	FGF21, Fibroblast growth factor 21, FGF-21
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	Q9NSA1
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Clone ID	Iv0067
Note	For research use only. Not suitable for clinical or therapeutic use.