

# Human IL10, InVivo Block/Neutralizing Antibody


## Overview

---

|                              |                                                                                                                                                              |
|------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No.</b>           | BGT-ANT-37359                                                                                                                                                |
| <b>Host species</b>          | Mouse                                                                                                                                                        |
| <b>Species reactivity</b>    | Human                                                                                                                                                        |
| <b>Form</b>                  | Liquid                                                                                                                                                       |
| <b>Storage buffer</b>        | 0.01M PBS, pH 7.4.                                                                                                                                           |
| <b>Purity</b>                | >95% as determined by SDS-PAGE.                                                                                                                              |
| <b>Clonality</b>             | Monoclonal                                                                                                                                                   |
| <b>Isotype</b>               | IgG1, kappa                                                                                                                                                  |
| <b>Applications</b>          | ELISA, Neutralization                                                                                                                                        |
| <b>Target</b>                | IL-10, IL10, Cytokine synthesis inhibitory factor, CSIF, Interleukin-10                                                                                      |
| <b>Purification</b>          | Protein A/G purified from cell culture supernatant.                                                                                                          |
| <b>Endotoxin level</b>       | Please contact with the lab for this information.                                                                                                            |
| <b>Accession</b>             | P22301                                                                                                                                                       |
| <b>Stability and Storage</b> | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |
| <b>Clone ID</b>              | Iv0025                                                                                                                                                       |
| <b>Note</b>                  | For research use only. Not suitable for clinical or therapeutic use.                                                                                         |

## Images

---



SEC-HPLC detection for Human IL10, InVivo  
Block/Neutralizing Antibody

Bioactivity